



USER MANUAL



WE HAVE YOU COVERED

Our passionate team brings you innovative, quality mobility solutions. Our products are reliable, sleek, and designed with simplicity so that you can rest assured that they work when you need them.

Thanks to our out of the box thinking, we are vigilant to meet your mobility needs with our simple and practical solutions. We recognize the important role that safety plays in consumer choices and are proud to offer a line of products that maintain the safety features found in your vehicle. Backed by third party verification for compliance to applicable Federal Motor Vehicle Safety Standards and Canadian Motor Vehicle Safety Standards, our products are engineered for safety and durability.

Our commitment to quality and safety along with our genuine desire to meet your mobility needs makes for a product we are sure you will be proud to own.





HI-LIFT

USER GUIDE

To load your wheelchair

The **HI-LIFT** is an optional wheelchair lift designed to lift manual rigid or folding wheelchair into your minivan via the sliding door.

Follow the next few steps:

- Once you have transferred out of your wheelchair, open the sliding door of your minivan to gain access to your **HI-LIFT**. Retrieve your hand held pendant.
- Press and hold the **OUT** button until the HI-LIFT has completely pivoted out of your vehicle.
- Press and hold the **DOWN** button until the strap of the **HI-LIFT** has lowered enough to reach the seat of your wheelchair. Remove your cushion and hook the strap to the loop that has been installed under the cushion of your wheelchair.



USER GUIDE (CONTINUED)

To load your wheelchair

- Press and hold the **UP** button until your wheelchair is high enough to enter the vehicle. When you release the button, the HI-LIFT will stop moving. Press and hold the **IN** button until your wheelchair is completely inside the vehicle. For this motion to go smoothly, we strongly suggest you keep a hand on your wheelchair (as pictured below) to help guide the chair into place. This will also prevent your wheelchair from sliding from side to side. Press and hold the **DOWN** button until your wheelchair touches the floor inside your vehicle. We suggest you then fasten your wheelchair in place.
- Stow the hand held pendant and close the vehicle's sliding door.

WARNING

WHEN OPERATING THE **HI-LIFT**, ALWAYS KEEP A HAND ON THE WHEELCHAIR TO PREVENT IT FROM SWINGING FROM SIDE TO SIDE.



⚠ WARNING: This product can expose you to chemicals including benzene, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.

USER GUIDE (CONTINUED)

To unload your wheelchair

Follow the next few steps:

- Open the sliding door and retrieve the hand held pendant.
- Once you have unfastened your wheelchair, press and hold the **UP** button until your wheelchair has lifted off the floor of the vehicle enough to clear the doorsill. Press and hold the **OUT** button until your wheelchair has come completely out of your vehicle. For this motion to go smoothly, we strongly suggest you keep a hand on your wheelchair (as pictured below) to help guide the chair out of the vehicle. This will also prevent your wheelchair from sliding from side to side. Press and hold the **DOWN** button until your wheelchairs has completely lowered to the ground.
- Unhook the **HI-LIFT** strap from your wheelchair. Press and hold the '**UP**' button until the strap of the **HI-LIFT** is fully stowed. Then, press and hold the **IN** button until the **HI-LIFT** has pivoted into the vehicle. Replace the pendant and close the door.





3 YEAR LIMITED WARRANTY

HI-LIFT

Adapt Solutions # 1 Ltd. (Adapt Solutions), warrants to the original purchaser of a **HI-LIFT** that the equipment is free from defects in material and workmanship for a period of three years from date of purchase.

During the first year of the warranty, Adapt Solutions will supply the replacement parts as well as a pre-set monetary amount (determined by Adapt Solutions) for the repair if a defect in material or workmanship is discovered. After the initial year of this warranty, only parts and components are covered. This warranty does not cover labour and other services after the initial year. Freight and other related repair charges will be the responsibility of the original purchaser.

The only remedy for a defect in one of Adapt Solutions products shall be the repair or the replacement, at the discretion of Adapt Solutions, of the defective part or component. If repair or replacement is not commercially practical or cannot be timely made, Adapt Solutions may decide to refund the purchase price of the equipment instead of repairing or replacing the original equipment.

In no event shall Adapt Solutions be responsible for indirect, incidental or consequential damages, whether such damages arise from claims based on contract, warranty, tort (including negligence), strict liability or product liability.

All implied warranties, including any warranty of merchantability or fitness for a particular purpose, are limited in their duration to the length of the warranty stated above for the affected component.

This warranty is to the original purchaser only, and excludes product damage due to installation error, product misuse, product abuse, accidents, physical damage, damage in shipment, modifications not made by Adapt Solutions, or repairs undertaken by anyone other than authorized distributors.

TO OBTAIN WARRANTY PARTS OR REIMBURSEMENT

You must follow these procedures:

1. Obtain warranty authorization by calling your local Adapt Solutions dealer or Adapt Solutions directly at 1-866-641-0419.
2. Return the faulty Adapt Solutions component/equipment, freight prepaid, to the address provided by your Adapt Solutions dealer or Adapt Solutions with proof of purchase indicating the date purchased.

Adapt Solutions will pay for shipping back to the purchaser within the continental United States and Canada if a defect in materiel or workmanship is discovered. Return freight and repair charges will be the responsibility of the purchaser if the problem is not covered by warranty.

This warranty gives you specific legal rights, and you may also have other rights that vary from state to state/ province to province. Adapt Solutions specifically does not authorize any person to extend the time or scope of this warranty.

For further information regarding this limited warranty, please contact us by calling 1-866-641-0419 or by writing at the following address:

**Adapt Solutions / Warranty Department
145 Damase-Breton
St-Lambert-de-Lauzon, Québec
G0S 2W0 Canada**

WARRANTY REGISTRATION

Adapt Solutions #1 Ltd located in St-Lambert-de-Lauzon, Québec, Canada, warrants this product to be free from any manufacturer's defect and workmanship for a period of three years from the date of purchase. This product has to be installed by an authorized dealer and has to be operated properly to be warranted. This warranty is limited to the original purchaser and does not cover defects on the vehicle. The warranty registration must be completed and returned within 30 days of purchase. This warranty does not cover abnormal wear or damages resulting from an accident, misuse, neglect or if the product has been altered in any way or aspect. Defective parts must be returned by an authorized dealer, who will have received a prior return authorization from Adapt Solutions #1 Ltd. All warranted parts shall be sent to Adapt Solutions #1 Ltd prepaid. Parts will be replaced or repaired at Adapt Solutions #1 Ltd.

Date installed: _____

Dealer: _____

Dealer signature: _____

TO KEEP FOR YOUR RECORDS

(Place serial number here)

WE ARE THE SOLUTION

Please contact your dealer for immediate assistance

Dealer:

TF 1.866.641.0419 | tech@adaptsolutions.com



145 Damase-Breton | Saint-Lambert, Qc G0S 2W0

TF 1.866.641.0419 | **Fax** 418.889.9838 | **P** 418.889.0419

www.adaptsolutions.com

